

SAMPLE PREPARATION

FL-0152

FLEC Carpet Test Plate, 25 mm



Overview

The FL-0152 Carpet Test Plate holds carpet and textile samples up to 25 mm thick for emission measurements with the FLEC Cell. FL-0151 adapters can be used in 1 mm increments to keep thinner sample surfaces level with the base of the cell. FLEC Cell is not included.

Technical overview

Product code	FL-0152
Sample thickness	Up to 25 mm
Samples	Carpets and textiles
Height adjustment	FL-0151 adapters in 1 mm increments
FLEC Cell	Not included

Related equipment

FL-0001 FLEC Cell

FL-0151 FLEC Test Plate Adapters, 5 pcs

[Product page and quotation enquiry](#)