


CORE SYSTEMS

FL-2002

FLEC Sub-Unit, 152 mm ID



Chamber base lowered



Chamber base raised



FL-2001 shown with FLEC Cell

The FL-2002 FLEC Sub-Unit supports emission sampling from non-planar, porous, liquid or sticky materials. Its 152 mm internal diameter accommodates a standard Petri dish. The internal height can be adjusted from 0 to 100 mm for dynamic or static headspace measurements, and the removable legs allow use on a heating plate or in an oven when required by the test procedure.

Specifications

Product code	FL-2002
Weight	13.5 kg
Dimensions	Height: 220 mm; outside diameter (OD): 200 mm
Chamber diameter	ID = 152 mm
Chamber height	0-100 mm
Sample container	Accommodates a standard Petri dish
Samples	Non-planar, porous, liquid or sticky materials
Mounting	Removable legs; heating plate or oven use when required by the procedure

[Product page and quotation enquiry](#)